



Our Scalars

Compare our scalars and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 18 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours.





- High-polish mirror finish tips.
- Manufactured within the guidelines of the FDA quality assurance system for good manufacturing practices.
- Highest grade German stainless steel, razor sharp, durable tips. Rockwell hardened to 60 RHc.
- 27 of the most popular designs from graceys to McCalls to columbias, and many european patterns.
- All instruments are guaranteed against defects in material and workmanship.

- All instruments bear the ce mark, emblematic of the european community's highest standards.
- Customer satisfaction guarantee, too!
- Less than 20 grams ... Lightweight for less hand fatigue.
- Manufactured to the most exacting standards so the 50th instrument you buy is exactly like your first.
- You decide ... Why pay more for the same high-quality instruments?

HANDLES FOR MOUTH MIRRORS

- Ручки для зеркал
- Manches pour miroirs à bouche
- Mundspiegelgriffe

- Ручки для дзеркал
- Mangos para espejos de boca
- مقابض لمرآيا الفم



DE-377 Hollow Mirror Handle Stem Extra Light (7.5mm) Length 5"



DE-378 Mirror Handle Stem Octagonal (6mm) Length 5"



DE-379 Mirror Handle Stem Round (6mm) Length 5"



DE1-380 Solid 12cm



DE2-380 Solid



DE3-380 Mirror Handle Stem with Rule 35mm (6mm) Length 5"



DE-380 Mirror Handle Stem with Rule 35mm (6mm) Length 5"



- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten

- Пинцети
- Pinzas para algodón
- ملاقط



DE-381
Surgical
11.5cm



DE-382
London College
15cm



DE-383
Semken Taylor
Str.



DE-384
Semken Taylor
Cvd.



DE-385
Meriam
16cm



DE-386
Dressing
14cm



DE-387
Semken Taylor
1x2 Straight



DE-388
Semken Taylor
1x2 Curved



DE-389
Gilles
15cm 1x2



DE-390
14cm 1x2



DE-391
Perry
13cm



DE-392
Atraumatic
15cm Cvd



DE-393
Atraumatic
15cm Str.



DE-394
College with Lock
15cm Cvd.



DE-395
Gerald
17.5cm



DE-396
with lock
15cm



DE-397
Atraumatic
17.5cm

TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten

- Пинцети
- Pinzas para algodón
- ملاقط



DE-398
Adson-TC
12cm



DE-399
Adson-TC
12cm 1x2



DE-400
Adson-TC
15cm 1x2



DE-401
Adson
12cm



DE-402
Adson
12cm 1x2



DE-403
Adson
15cm 1x2



DE-404
Adson Brown
12cm



DE-405
Crane Kaplan
15cm



DE-406
Crane Kaplan
15cm



DE-407
Crane Kaplan
15cm



DE-408
Hemingway
15cm



DE-409
Miller
15cm Str.



DE-410
Miller
15cm Cvd.



DE-411
Dressing
Tweezer
14.5cm



DE-412
College
Curved
Tip 15cm



DE-413
Maryam
Tweezers
15cm



DE-414
College
Tweezer
15cm



DE-415
Perry
Tweezer
12.5cm



DE-416
Self-locking
Tweezer
15cm

- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين

Single Ended Probe Solid Handle 4.5mm Length 150mm




Fig. 2
Normal Handle

Single Ended Probe Hollow Handle 8.5mm Length 150mm




Fig. 2
Anatomic Handle


EXPLORERS

- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين

Probes Single End Solid Handle 6mm, length 150mm



Probes Single Ended 6mm, length 150mm



Examination kit Set of 3 Pieces



Examination Set of 4 Pieces




Examination Set of 4 Pieces



**EXPLORERS INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM**


- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين




EXPLORERS

- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين

Pointed Probes Single Ended Handle 6mm, length 150mm



Pointed Probes Double Ended Handle 6mm, length 170mm



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie

- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة

PERIODONTAL SET OF 10 PIECES SOLID HANDLE 6MM LENGTH 170MM




EXPLORERS

- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين

WILLIAM PROBE SET DOUBLE ENDED WITH HOLLOW HANDLE 9.5MM, LENGTH 170MM
CERAMIC COATED HANDLE


Proper care must be given during cleaning, sharpening and sterilization to preserve design and integrity and to prolong usefulness. Do not expose instruments to sterilization temperatures above 350 F/177 C or to higher than recommended chemical concentrations.



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie

- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة

PERIODONTAL INSTRUMENTS SET OF 11 PCS SOLID HANDLE 6MM, LENGTH 170MM



CURETTES

- Кюрети
- Curetas
- Küretten

- Кюрети
- Curette
- كرتس

GRACEY SET OF 7 PCS HOLLOW HANDLE 9.5MM, LENGTH 170MM



• Кюрети
• Curetas
• Küretten

• Кюрети
• Curette
• كرتس

GRACEY SET OF 7 PCS
HOLLOW HANDLE 9.5MM, LENGTH 170MM WITH SILICON GRIP



DE-524
GRACEY
1/2



DE-525
GRACEY
3/4



DE-526
GRACEY
5/6



DE-527
GRACEY
7/8



DE-528
GRACEY
9/10



DE-529
GRACEY
11/12



DE-530
GRACEY
13/14



CURETTES

- Кюрети
- Curetas
- Küretten

- Кюрети
- Curette
- كرتس

GRACEY MINI FIVE


HOLLOW HANDLE 9.5MM, LENGTH 170MM
CERAMIC COATED HANDLE



DE-531 Gracey 1/2
DE-532 Gracey 3/4
DE-533 Gracey 5/6
DE-534 Gracey 7/8
DE-535 Gracey 11/12
DE-536 Gracey 13/14

GRACEY AFTER FIVE

HOLLOW HANDLE 9.5MM, LENGTH 170MM
CERAMIC COATED HANDLE



DE-537 Gracey 1/2 AF
DE-538 Gracey 3/4 AF
DE-539 Gracey 5/6 AF
DE-540 Gracey 7/8 AF
DE-541 Gracey 11/12 AF
DE-542 Gracey 13/14 AF

Proper care must be given during cleaning, sharpening and sterilization to preserve design and integrity and to prolong usefulness. Do not expose instruments to sterilization temperatures above 350 F/177 C or to higher than recommended chemical concentrations.



- Кюрети
- Curetas
- Küretten
- Кюрети
- Curette
- كرتس

**GRACEY CURETTE
HOLLOW HANDLE 10MM, LENGTH 170MM**

Curetters Set of 8 Pieces



FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء التصكوك

**FILLING INSTRUMENTS (Gold Pointed)
SOLID HANDLE 6MM, LENGTH 170MM**




EXPLORERS INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM

- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين

PLUGGERS

2/4
2/3
2/2
2/1
2/0
1/9
1/8
1/7
1/6
1/5
1/4
1/3
1/2
1/1
1/0
0/9
0/8
0/7
0/6
0/5
0/4
0/3
0/2
0/1cm
0
CENTIMETERS




PERIODONTAL PROBES

- Зонды пародонтологические
- Sondes parodontales
- sonda parodontali
- зонди пародонтологічні
- Sondas parodontales
- المسابر اللثة

SPREADERS


HEAT CARRIERS



- Зонды пародонтологические
- Sondes parodontales
- sonda parodontali
- зонди пародонтологічні
- Sondas parodontales
- المسابر اللثة

ENDODONTICS INSTRUMENTS

ANTERIOR



A QUALITY ORIENTED COMPANY
Dentech
Surgical (Pvt) Ltd.

LUKS



DE-587
8A
29mm

PERIODONTAL PROBES


- Зонды пародонтологические
- Sondes parodontales
- sonda parodontali
- зонди пародонтологічні
- Sondas parodontales
- المعابير اللثة

ENDODONTICS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

POSTERIOR



RETRO FILLING PLUGGERS




- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)


- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

POSTERIOR



UNIVERSAL FILLING INSTRUMENTS



FILLING INSTRUMENTS


- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BLACK

PLUGGERS/CONDENSERS

Hollenback




- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)


- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

PLUGGERS/CONDENSERS

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM



PLUGGERS/CONDENSERS




FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

Tanner Filling Instruments



DE-649
2T
Tanner

DE-650
3T
Tanner



DE-651
9/10
Back Action



DE-652
11/12
Back Action




DE-653
153
English Pattern
Serrated



DE-654
153
English Pattern
Serrated

BURNISHERS

BALL



DE-655
18
#41, #6

DE-656
26/27S
#41, #6

DE-657
27/29
#41, #6, #7, #8

BEAVERTAIL



DE-658
2

DE-659
2
Double-End

ACORN



DE-660
21
#41, #6

DE-661
21B
#41, #6, #8

LADMORE



DE-662
2
Ladmore

DE-648
1T
Tanner



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BURNISHERS

PERIODONTAL FILES




FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE/PLASTIC FILLING INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM

BURNISHER FILLING INSTRUMENTS




UNIVERSAL FILLING INSTRUMENTS



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM



AMALGAM CARVERS

- Инструменты для реставрации
- Instruments pour amalgame
- Amalgampistolen
- Инструменти для реставрації
- Instrumentos para amalgama
- تقطيع الملعف

AMALGAM INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM



- Инструменты для реставрации
- Instruments pour amalgame
- Amalgampistolen

- Инструменти для реставрації
- instrumentos para amalgama
- تقطیع المنظم

AMALGAM INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

Excavators Instruments Set of 11 Pieces



PERIODONTAL INSTRUMENTS

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- Instrumentos para periodoncia
- المصكوك اللثة

SICKLE SCALERS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

SICKLE SCALERS



- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie

- Інструменти пародонтологічні
- Instrumentos para periodoncia
- الصكوك اللثة

PERIODONTAL INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

JACQUETTE SCALERS

CRANE-KAPLAN



PERIODONTAL INSTRUMENTS

- Инструменты пародонтологические
- Instruments pour la parodontologie
- Instrumente für die Parodontologie
- Інструменти пародонтологічні
- Intrumentos para periodoncia
- الصكوك اللثة

PERIODONTAL INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM SCALER/CURETTE SETS



- Кюрети
- Curetas
- Kuretten
- Кюрети
- Curette
- کرتس

SCALER/CURETTE SETS

GOLDMAN-FOX SET




CURETTES

- Кюрети
- Curetas
- Küretten

- Кюрети
- Curette
- كرتس

GRACEY CURETTE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

GRACEY CURETTES



DE-757
9/10†
(shown)
#2, #4, #6
#7, #8, #9

DE-758
9/10
Rigid
#2, #4, #6
#7, #8

DE-759
11/12†
(shown)
#2, #4, #6
#7, #8, #9

DE-760 **DE-761**
11/12 11/12
Rigid Extra Rigid
#2, #4, #6 #4, #6 #7,
#7, #8, #9 #8

DE-792
13/14†
(shown)
#2, #4, #6
#7, #8, #9

DE-793 **DE-764**
13/14 13/14
Rigid Extra Rigid
#2, #4, #6 #4, #6 #7,
#7, #8, #9 #8

DE-765
15/16†
(shown)
#2, #4, #6
#7, #8, #9

DE-766
15/16
Rigid
#2, #4, #6
#7, #8, #9

DE-767
17/18
#2, #4, #6
#7, #8, #9

DE-768
17/18
Rigid
#2, #4, #6
#7, #8, #9

- Кюрети
- Curetas
- Küretten
- Кюрети
- Curette
- كرتس

**CURETTE INSTRUMENTS
HOLLOW HANDLE 8.5MM, LENGTH 170MM**



DE-769
Hemingway
2mm, 2mm

DE-770
Hemingway
3mm, 3mm

DE-771
Hemingway
4mm, 4mm

**CURETTE INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM**



DE-772
Hemingway
2mm, 2mm

DE-773
Hemingway
3mm, 3mm

DE-774
Hemingway
4mm, 4mm

**CURETTE INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM**



DE-775
Hemingway
2mm, 2mm

DE-776
Hemingway
3mm, 3mm

DE-777
Hemingway
4mm, 4mm

EXCAVATORS

- Экскаваторы
- Excavateurs
- Excavatoren
- عسكافطوري
- Excavadores
- حفارات

EXCAVATORS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

SPOONS



- Экскаваторы
- Excavateurs
- Excavatoren

- Экскаватори
- Excavadores
- حفارات

OVAL SPOONS/DARBY-PERRY

EXCAVATORS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

EXCAVATORS



EXCAVATORS

- Эскаваторы
- Excavateurs
- Excavatoren
- عسكافاتوري
- Excavadores
- حفارات

EXCAVATORS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BLADES




• Экскаваторы
• Excavateurs
• Excavatoren

• Екскаватори
• Excavadores
• حفارات

EXCAVATORS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

ENGLISH PATTERNS



DE-802
Goldstein
Interproximal
Excavator

DE-803
243

DE-804
244

DE-805
245

DE-806
246

DE-807
125/126
#41, #6

DE-808
127/128
#41, #6

DE-809
129/130
#41, #6

DE-810
131/132
#41, #6


DE-811
153/154

FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

HATCHETS

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

CARVERS

RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM



Tip Size: 1/1

Solid Handle Size 4.25"

DE-827
Inter-proximal
Off-Angle

DE-828
1/2-3
Hollenback
#41, #6,
#7, #8

DE-829
3S
Hollenback
#41, #6, #7

DE-830
3
Hollenback

DE-831
8
Wiland
Minn.IPC1
#41, #6, #7

TANNER

WALL

WARD

DE-826
Inter-proximal
#41, #6, #7

DE-832
3
Tanner

DE-833
4
Tanner

DE-834
5
Tanner

DE-835
6
Tanner

DE-836
1
Wall

DE-837
3
Wall

DE-838
1
Ward

DE-839
2
Ward

AMALGAM CARVERS

- Инструменты для реставрации
- Instruments pour amalgame
- Amalgampistolen
- Инструменти для реставрації
- Instrumentos para amalgama
- تقطيع الملعف


RESTORATIVE INSTRUMENTS

ANDREW

SOLID HANDLE 6MM, LENGTH 170MM


CLAPPISON/FRAHAM

HOLLENBACK



DE-841
2
Andrew

DE-842
3
Andrew



DE-843
4
Clappison

DE-844
2/3
Frahm
#41, #6

DE-845
104



DE-846
1 Hollenback
2 Modified
Cleoid Small

DE-847
1 Hollenback
2 Modified
Cleoid Large


HOLLENBACK

THARP



DE-848
1
Hollenback


DE-849
2
Hollenback



DE-850
L7
Levy

DE-851
SHO-A
Shoshan

DE-852
18
Sprengel



DE-853
1
Tharp

DE-854
2
Tharp

DE-840
1
Andrew

- Инструменты для реставрации
- Instruments pour amalgame
- Amalgampistolen

- Инструменты для реставрації
- instrumentos para amalgama
- تقطیع المنعم


RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

CARVERS

AMALGAM FILES



HOES

Forming Line Angles on Anterior Preparations.



DISCOID-CLEOIDS

DISCOID-CLEOIDS




FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

RESTORATIVE INSTRUMENTS
SOLID HANDLE 6MM, LENGTH 170MM

MARGIN TRIMMERS




- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente

- Інструменти для реставрації
- Instrumentos para modelar
- أدوات تلطيع

P.K. THOMS WAXING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

P.K. THOMAS WAXING INSTRUMENTS



DE1-884
P.k. Thomas Waxing
Instruments
Kit



COMPOSITE FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

Composite Gold Tip Instruments



DE1-892

DE2-892

DE3-892

DE4-892

DE5-892

DE6-892

DE7-892

DE8-892

DE9-892

COMPOSITE FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك


Composite Gold Tip Instruments



FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC



DE21-892

DE22-892

DE23-892

DE24-892

DE25-892

DE26-892



DE27-892
Composite Instruments Blue
Titanium with Cassettes

- Прочие товары
- Divers produits et services
- Sonstige Produkte

- Інші товари
- Varios Productos
- المنتجات المتنوعة



FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

CHISELS

HOES



CEMENT SPATULAS

- Шпатели для цемента
- Spatules à ciment
- Zementspatel

- Шпателі для цементу
- Espátulas para cemento
- ملاعق الأسمنت



DE-893
00



DE-894
0



DE-895
1



DE-896
2



DE-897
3



DE-898
4



DE-899
3
Weston



DE-900
1
medium cement



DE-901
1



DE-902
Le Cron



DE-903
Le Cron



DE-904
Le Cron



DE-905
MSY



DE-906
MSY



DE-907
Le Cron



DE-908
Hylin



DE-909
Hylin



DE-910
Hylin



DE-911
Zahle



DE-912
Zahle



DE-913
Zahle



DE-914
Gritman

WAX KNIVES

- Ножи для воска
- Couteaux a cire
- Wachsmesser
- Ножі для воску
- Cuchillos para cera
- السكاكين الشمع

ALGINATE, PLASTER SPATULAS, WAX KNIVES PLASTER FOR KNIVES & PLASTER SHEAR



DE-915
Wooden Handle
Fig. 5
21.5 cm



DE-916
Wooden Handle
Fig. 7
19.5 cm curved



DE-917
Wooden Handle
Fig. 8
20 cm curved



DE-918
Wooden Handle
Fig. 9
22.5 cm curved



DE-919
Hook Knives



DE-920
Fig. 1
Lessmann
12.5 cm



DE-921
Fig. 2
Lessmann
17 cm



DE-922
Fig. 1
with small spoon
12.5 cm



DE-923
Fig. 2
with large spoon
17 cm



DE-924
Bacolite Handle
Wax Knives,
Small



DE-925
Bacolite Handle
Wax Knives,
Large



DE-926
Wooden Handle
Wax Knives,
Small



DE-927
Wooden Handle
Wax Knives,
Large



DE-928
Plaster Knife



DE-929
Fig. 1
17 cm



DE-930
Fig. 3
16 cm



DE-931
Fig. 4
16 cm



DE-932
Special Inox
19.5 cm

- Артикуляторы
- Artocilar
- Artikulations

- Артикулятори
- Articuler
- ارنیکو لائوس



DE-933

Technical Datas:

Intercondylar distance	120 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	10°
Distance between mounting plates	85 mm.
Weight	400 gr.
Opening (item 5000)	125°
Opening (item 5010)	180°



DE-934

Technical Datas:

Intercondylar distance	110 mm.
condylar inclination	fixed at 30°
Bennet angle	fixed at 7°
Incisal table inclination	flat
Distance between mounting plates	85 mm.
Weight	450 gr.
Opening	180°



DE-935

Stainless Steel Finishing

DE-936

Royal Blue



DE-937



DE-938



DE-939



DE-940



DE-941



DE-942



DE-943



DE-944



DE-945



DE-946
Nakawalni
Dental Laboratory
Instrument